



FALCON™ BATTERY MINDER™

CLIENT PROBLEM

Clients seem unable to remember that wheelchair and scooter batteries need to be recharged REGULARLY – literally every day without waiting for the batteries to wear down. Clients get wheelchair batteries confused with their electronic device batteries which work totally differently. When they don't charge regularly, two things happen and both are bad. The first is that performance of the chair suffers and the second is that the battery dies a premature death – leading to the necessity of purchasing more batteries than you should.

THE FALCON™ SOLUTION

The Falcon™ Battery Minder™ attaches directly to the batteries and monitors charge level. When the battery reaches 30% of charge remaining, the alarm repeats a programmable alert every 30 seconds reminding the client to recharge the batteries. If this alert is ignored then when 20% of the charge remains, a more strident alert issues every 30 seconds, until the charger is engaged. Since the alarm message is easily recorded, the clinician can determine the appropriate message including the language of choice.



SPECIFICATIONS

- Hard wired into the 24V system
- Repeats message every 30 seconds
- Clinician can record appropriate messages.

ORDERING INFORMATION

Part number: BMAV
List Price: \$595 each