



## Patient Alarms Price List - It's as easy as A B C

Sales Order No.		
Mark for	Buyer	Purchase Order
Product Code:	Phone No. ( )	Fax No. ( )
<b>BILL TO:</b>		<b>SHIP TO:</b>
Company	Company	
Attention	Attention	
Address	Address	
City/State/Zip	City/State/Zip	
Phone no. ( )	Phone no. ( )	

### A BATTERY MINDER FEATURES

THE FALCON™ BATTERY MINDER™ ATTACHES DIRECTLY TO THE BATTERIES AND MONITORS CHARGE LEVEL. WHEN THE BATTERY REACHES 30% OF CHARGE REMAINING, THE ALARM REPEATS A PROGRAMMABLE ALERT EVERY 30 SECONDS REMINDING THE CLIENT TO RECHARGE THE BATTERIES. IF THIS ALERT IS IGNORED THEN WHEN 20% OF THE CHARGE REMAINS, A MORE STRIDENT ALERT ISSUES EVERY 30 SECONDS, UNTIL THE CHARGER IS ENGAGED. SINCE THE ALARM MESSAGE IS EASILY RECORDED, THE CLINICIAN CAN DETERMINE THE APPROPRIATE MESSAGE INCLUDING THE LANGUAGE OF CHOICE.

#### Specifications

- Hard wired into the 24V system
- Repeats message every 30 seconds
- Clinician can record appropriate messages.



<input type="checkbox"/>	BMAV	Falcon Rehab Battery Minder	Qty:	@ \$595 each	Total	\$
--------------------------	------	-----------------------------	------	--------------	-------	----

### B SHIFT MINDER FEATURES

FALCON™ DEVELOPED THE SHIFTMINDER™ ALARM AND DATA COLLECTION UNIT THAT CAN HELP TRAIN NEW SYSTEM USERS, PROVIDE A SAFEGUARD FOR CLIENT'S WITH IMPAIRED COGNITIVE FUNCTION WHO MAY NOT ALWAYS REMEMBER TO MAINTAIN A PRESCRIBED SHIFT REGIMEN, AND PROVIDE DOCUMENTATION OF CLIENT BEHAVIOR FOR CLINICAL REVIEW.

#### Specifications

- Runs on 9-Volt Battery
- Data is exportable to MS Excel
- Can hold 2 months of client activity
- Available with Tone or Voice Alarm
- Maintains alarm function for over 670 hours before battery needs replacement



<input type="checkbox"/>	SMAT	Falcon Rehab Shift Minder – Tone Alarm	Qty:	@ \$595 each	Total	\$
<input type="checkbox"/>	SMAV	Falcon Rehab Shift Minder – Voice Alarm	Qty:	@ \$1095 each	Total	\$

### C CUSHION MINDER FEATURES

THE CUSHIONMINDER™ ALARM IS A WAY TO ENSURE THE CORRECT POSITION OF A CUSHION IN A CLIENTS' WHEELCHAIR (NOT UPSIDE DOWN OR BACKWARDS) SINCE THE ALARM MESSAGE IS EASILY RECORDED, THE CLINICIAN CAN DETERMINE THE APPROPRIATE MESSAGE INCLUDING THE LANGUAGE OF CHOICE.

#### Specifications

- Hard wired into the 24V system
- Repeats message every 30 seconds
- Clinician can record appropriate messages.

<input type="checkbox"/>	CMAV	Falcon Rehab Cushion Minder	Qty:	@ \$595 each	Total	\$
--------------------------	------	-----------------------------	------	--------------	-------	----

Falcon Rehabilitation is not responsible for typographical errors. Packaging, pricing, warranties and products are subject to change without notice.