



Falcon wanted to give you a simple reminder of all of the "BEAR"iatric products we now make.

WE made it cute so you will want to save it.



Manual Chairs		Power Seating	
MTC	Photon	Tilt	
<p>The manual wheelchair series includes tilt, recline, tilt and recline with manual or power functions.</p> <p>These products can support up to 600 pounds and are available in innumerable widths, depths and back heights.</p>	<p>The Photon offers a light weight chair for the larger, active user. The Photon can accommodate up to 300# and can be ordered in widths and depths of 20 x 20. An 18x18 Photon w/ pneumatic wheels weighs 19.75 #.</p>	<p>The power seating is available in tilt, recline, tilt and recline with all the standard power options. Weights and configurations up to 650 pounds, and widths of 28" wide and 22" deep are considered standard.</p> <p>Outside these parameters the product is considered custom.</p>	
Power Accessories		Non - Power Accessories	
Scissors Lift - SKU	Windswept legrest	Offset Hanger	Insert Back Posts
<p>The scissor seat lift is designed to work with most bases and most seats currently on the market. The 4 points of stability provide a secure seat elevation of up to 10". It works equally well with rehab seating or the captain's seat used in standard power mobility.</p>	<p>Heavy duty multi axis legrests are ideal for the client whose stomach distends and pushes the legs apart necessitating a windswept legrest to follow the contours on the leg outside the parallel of the wheelchair's profile. You can open up the hip angle in 6 second degrees of arc to accommodate the client's needs.</p>	<p>Offset legrest hangers is a fixed approach for the windswept client. It provides additional space between the legrests, and allows the use of other manufacturer's legrests.</p>	<p>Inset, backrests are highly recommended for the bariatric client who carries much of the weight in and below the hips. It allows you to space the arms closer together thus reducing the overall width of the wheelchair. Anytime you are confronted with a seat width of 20" or more it should be considered.</p>

Specialty Items			
Power Panniculift	Lateral Support	Bariatric Arm	Prone Carts
			
<p>Power Panniculift fits into the "who knew" category. For patients who'd lost large amounts of weight, but were still in process, the distended tissue prevents sitting upright, interferes with breathing and digestion, and makes ADL's difficult or impossible. The power Panniculift is a solution for the user waiting to reach ideal weight and the panniculectomy that will finally correct the problem.</p>	<p>The bariatric requirements for lateral support include its ability to support the side weight if the client leans significantly left or right. The Falcon lateral supports will hold up to 350# of side pressure. A product not made out of steel and not having positive teeth to hold the position won't keep the client in proper alignment. We err on the side of too strong.</p>	<p>Our cantilevered arms can support the full weight of the client's push. The arms are designed to help when the hip extends beyond the width of the chair because there is no post needed to support the arm, This helps keep doorways navigable.</p> <p>These arms are compatible with most systems currently on the market.</p>	<p>Most of our family of Prone Carts can be made to accommodate the rare client who has both bariatric considerations as well as being in a long-term prone position.</p> <p>Prone carts can be made in manual or power assisted configurations, in a variety of lengths and with additional options. For more specifics please see our Prone Cart Product Information Guide.</p>



A Note from Mr. Moose: "Here at Falcon we look at the fairy tale of Jack Spratt and his bariatric wife as more of a discretionary tale. Remember, it doesn't take that many fast food hamburgers to start sliding down that precipitous slope".